



# Ultima 5

DIGITAL TENS STIMULATOR

## Features

- ▶ *Dual Channel (4-lead) Digital TENS*
- ▶ *Two Selectable Waveforms: (Symmetrical biphasic or Monophasic)*
- ▶ *Five Treatment Modes (Constant, Burst I, Burst II, Modulated Width, Modulated Rate & Width)*
- ▶ *Large LCD display showing all parameters*
- ▶ *Last Setting Recall & Settings Lock*
- ▶ *Treatment Timer (15, 30, 45, 60, 90 minutes, or Continuous)*
- ▶ *Removable plastic belt clip*
- ▶ *1-year warranty*

## Technical Specifications

<b>Channel:</b> Dual
<b>Amplitude Maximum:</b> 150mA (peak value) across 500 ohm load
<b>Output Voltage:</b> 0-73 V across 500 ohm load
<b>Waveform:</b> Symmetrical Biphasic, or Monophasic
<b>Pulse Rate:</b> Adjustable, 1 - 150 Hz
<b>Pulse Width:</b> Adjustable, 50 - 250µS in steps of 10 µS Treatment
<b>Timer:</b> 15, 30, 45, 60, or 90 minutes, or Continuous
<b>Treatment Modes:</b> Constant, Burst I, Burst II, Modulated Width, Modulated Rate & Width
<b>Power Supply:</b> 2 AA alkaline batteries 1.5 V
<b>Dimensions:</b> 10 x 6.5 x 2.8 cm
<b>Weight:</b> 127.5 grams (including batteries)

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA  
1-800-363-STIM | Fax: 1-651-565-2410 | www.medi-stim.com

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43 pin lead wires, 2 each AA alkaline batteries, carrying case, and instruction manual

